# GAMINIO NAUDOJIMO INSTRUKCIJA USER'S MANUAL

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Applicable for the treatment of epicondylitis. Provides partially stabilization of tendons on the outside of the forearm (tendons attaching to external humerus). Especially recommended for tennis player and physiotherapists.

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**Composition:** 75% elastic polyurethane, 12.5% polyamide, 12.5% viscose.

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**Size.** In order the product performs its functions effectively correct size has to be chosen (see Table below). Orthosis should not slip down, strangle or cause any other unfavourable sensations.

	XS	S	M	L
cm	23-25	25-27	27-30	30-33

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### Directions for use:

1. Put the setting on the hand below the elbow joint in such a position that a clasp would be outside the hand, while the fixing pad would stress the muscle 2-3 cm. below the elbow flexure or below the outside crank web papilla.

2. Hook up the sticky strip, regulating the upper girth in such a position that a hand couldn't go up to sleep (See Picture 2).

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### Recommendations:

Consult your doctor on the expediency of wearing the orthosis and the recommended time of wearing it. Do not wear the orthosis if you notice any skin sensitisations.

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**The warranty** is provided in accordance with the procedure established by the laws. The warranty obligations shall be observed provided that during the said period the product is used for its direct purpose and its construction is not modified.

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**Thank you for choosing our product. We trust you'll feel better.**

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